

# Teeth Whitening

## Fact Sheet

### How does Rejuvenate! use light to whiten teeth?

The Rejuvenate! teeth whitening system starts with a patented whitening gel, applied directly to the surface of the teeth. The powerful whitening accelerator then produces an intense beam of pure LED light filtered through over 12,000 optical fibers and 30 layers of coated optical glass, illuminating the teeth at a wavelength proven to be the optimal range for tooth whitening. This light activates the whitening gel and rapidly accelerates the whitening process. The results are much more effective than home whitening products, without the sensitivity or discomfort, and the results are immediate.

### How long does the treatment take?

In less than an hour, tooth discoloration from years of staining foods, coffee, tobacco, cola and red wine is eliminated, leaving a beautiful, pearly-white smile.

### How much of an improvement will I see?

Average results in shade improvement are between 5 and 14 shades on the whitening scale.

### Is it safe?

Yes. Any potentially-harmful ultraviolet (UV) and infrared light is completely removed by the filtering system, leaving only a focused beam of cold light in the blue-cyan spectrum.

### Does it hurt?

The teeth whitening process is completely painless, and the results are immediate and long-lasting. A small percentage of clients may feel minor sensitivity in the first twelve hours following the treatment.

### How do I book treatment?

Call Rejuvenate! Med Spa today at 866.594.5476 for a free consultation. We will then arrange a treatment session that seamlessly fits into your schedule.



Unretouched clinical photos taken before and after treatment with the Rejuvenate! teeth whitening System. Results may vary.



**rejuvenate!**

med spa | cosmetic surgery | luxury spa

toll free: 866.594.5476

info@rejuvenate-medspa.com

rejuvenate-medspa.com